

# **Circumferential Correction of Severe Thoracolumbar Kyphosis by Utilizing an Articulating Rod and Lateral Interbody Device: A Technical Note**

Sumedh S. Shah, Malek Bashti, Manav Daftari, James Boddu, Nathaniel B. Dusseau II, Jason Liounakos and Timur Urakov

*Int J Spine Surg* published online 18 March 2025 https://www.ijssurgery.com/content/early/2025/03/17/8723

This information is current as of May 1, 2025.

**Email Alerts** Receive free email-alerts when new articles cite this article. Sign up at: http://ijssurgery.com/alerts



International Journal of Spine Surgery, Vol. 00, No. 0, 2025, pp. 1–5 https://doi.org/10.14444/8723 © International Society for the Advancement of Spine Surgery

# Circumferential Correction of Severe Thoracolumbar Kyphosis by Utilizing an Articulating Rod and Lateral Interbody Device: A Technical Note

# SUMEDH S. SHAH, MD<sup>1</sup>; MALEK BASHTI, BS<sup>1</sup>; MANAV DAFTARI, BS<sup>1</sup>; JAMES BODDU, MD<sup>1</sup>; NATHANIEL B. DUSSEAU, II, MS<sup>1,2</sup>; JASON LIOUNAKOS, MD<sup>1,3</sup>; AND TIMUR URAKOV, MD<sup>1</sup>

<sup>1</sup>Department of Neurological Surgery, University of Miami Miller School of Medicine, Miami, FL, USA; <sup>2</sup>College of Medicine and Life Sciences, The University of Toledo, Toledo, OH, USA; <sup>3</sup>Baptist Health Brain & Spine Care, Miami Neuroscience Institute, Baptist Health Medical Group, Miami, FL, USA

#### ABSTRACT

**Background:** Surgical management of severe thoracolumbar kyphosis (TLK) is challenging due to the necessity of achieving adequate sagittal realignment without increasing operative risk or compromising hardware integrity. Performing large corrective maneuvers during deformity correction may increase the risk of neurological complications, especially when correcting hyperkyphosis in the distal thoracolumbar spine. The present article describes the first use of a custom-made articulating rod to manipulate a patient's severe TLK into extension in combination with the expandable lateral interbody device to achieve circumferential deformity correction.

**Clinical Presentation:** A 38-year-old woman with a history of spinal trauma 10 years ago developed severe TLK ( $<70^{\circ}$  angulation) after failed posterior spinal stabilization. Due to her debilitating neurological symptoms and profound sagittal deformity, we performed a T9 to L2 navigation-assisted posterolateral fusion followed by a circumferential kyphotic correction utilizing a custom-made articulating rod (Globus Medical, Inc., Audubon, PA) with an expandable device placed in prone-lateral position at T11 to T12. The use of an articulating rod enabled us to safely manipulate the thoracolumbar spine into extension prior to permanent rod fixation. We were able to achieve approximately 40° of thoracolumbar correction. No immediate or late (at 2-year follow-up) postoperative medical or hardware-related complications were reported.

**Conclusion:** The use of an articulating rod with prone-lateral placement of an expandable interbody device allowed for a high degree of circumferential thoracolumbar deformity correction in a patient with severe post-traumatic TLK. Our results indicated the technical feasibility and success of utilizing this treatment strategy for high-grade TLK without medical- or hardware-related failure.

**Clinical Relevance:** The use of an articulating rod for correction of complex TLKs may benefit patient outcomes and improve safety.

Level of evidence: 4.

Novel Techniques & Technology

Keywords: articulating rod, kyphosis, lateral interbody fusion, prone-lateral, thoracolumbar deformity

### INTRODUCTION

Thoracolumbar kyphosis (TLK) represents a significant source of functional, neurological, and cosmetic morbidity.<sup>1,2</sup> Surgical management of severe TLK is challenging due to the necessity of achieving adequate realignment without increasing operative risk. However, in performing the necessary corrective maneuvers, the rate of neurological complication can be exceptionally high.<sup>3–5</sup> Etiology of neurological complications can range from blood flow and pressure aberrations, compression/traction of neural elements, or osteotomy-related injury.<sup>6</sup> Therefore, various techniques are employed to maximize corrective ability while decreasing operative risks, such as multistaged surgeries,

implementation of nontraditional approaches, and use of novel technologies for controlled spinal reduction to prevent neural injury.<sup>4</sup>

In the current article, we present a technical report of our patient undergoing circumferential correction of severe TLK ( $>70^\circ$ ) by combining a custom-made articulating rod to control compression/distraction forces on the neural elements and an expandable lateral interbody device. The articulating rod acts by capitalizing on a central pivot point that allows for dynamic maneuvers of the spinal alignment not possible with static rods and in situ benders. We believe that the collective use of such devices can allow for safe, durable, and significant thoracolumbar correction, particularly in highly sensitive spinal segments.

Downloaded from https://www.ijssurgery.com/ by guest on May 1, 2025





Figure 1. (A) Preoperative computerized tomography (CT) shows a severe thoracolumbar kyphosis due to wedge deformity causing posterior canal stenosis and a Cobb angle of 73° at the T11 to 12 apex. (B and C) Anteroposterior and lateral intraoperative fluoroscopy images showing posterolateral instrumented fusion with a left-sided approach lateral interbody device at T11 to T12 with integrated fixation. (D) Postoperative follow-up CT at 2 years with nearly 40° correction. No hardware issues were noted.

# CLINICAL PRESENTATION

A 38-year-old woman with a history of spinal trauma requiring previous lumbar fusion 10 years prior presented with complaints of severe positive sagittal deformity contributing to neurogenic claudication, back pain, and radiculopathy into the anterior thigh and knee bilaterally. Past medical history was significant for hypertension and obesity. While she was neurologically full strength and had no sensory deficits on physical examination, she exhibited a forward-leaning station and profoundly antalgic gait. Medications included gabapentin, methocarbamol, and buprenorphine. Conservative management including physical therapy, and chiropractic care was trialed without benefit. Computerized tomography and magnetic resonance imaging indicated a severe TLK with apex at T11 to T12 and a Cobb angle between T10 and L1 of 73° (Figure 1A). Evidence of prior failed fusion was noted with fracture and retained pedicle screw fragments at L2. Otherwise, available neuraxial imaging was benign.

#### **Operative Plan**

Given the degree of symptomatic kyphotic deformity and failure of conservative therapy, surgical intervention was recommended by way of T10 to L1 laminectomy, Ponte osteotomies at T11 to T12, placement of expandable lateral interbody device (Globus Medical Inc., PA, USA) at T11 to 12 in prone-lateral position, and T9 to L2 navigation-assisted posterior pedicle screw fixation. Considering the need to achieve significant sagittal deformity correction, we opted to use a custom-made articulating rod (Globus Medical Inc., PA, USA) to allow for controlled distraction and extension during



**Figure 2.** Schematic of an articulating rod using a commercially available 3-dimensional software used to create the device.

reduction maneuvers (Figure 2, Supplemental Video 1). The risks and benefits were discussed with the patient, and informed consent was obtained for operative management. The patient provided consent to use of their deidentified clinical information for publication. No institutional review board approval was necessary per institutional case report policy.

#### Surgical Management

Instrumentation and decompression with Schwab Grade 2 posterior osteotomies were performed. A standard superficial left-sided lateral approach was then utilized at the T11 to T12 level, tailored to the severity and anatomical presentation of the deformity. An anterior column release was performed in the lateral position using an anterior column blade to precisely transect the anterior longitudinal ligament, with surrounding structures protected throughout the procedure using a spec retractor. Due to the anatomical presentation of the deformity, the mobilization of the diaphragm and pleura was not required, and a retropleural approach was unnecessary. Segmental vessels were carefully ligated using bipolar cautery to prevent bleeding.

The articulating rod was then introduced and secured into the pedicle screws to sequentially expand the T11 to T12 disc space. This allowed for the placement of an expandable interbody device, inserted in the pronelateral position. Once adequate expansion was achieved, the articulating rod was adjusted into a slight extension, and the T12 to L2 screw caps were tightened. Reduction towers were placed on the top T9 to T11 screws to further correct the thoracolumbar deformity. No intraoperative neuromonitoring (IONM) changes were noted during these maneuvers.

After achieving the desired correction, the final rod was placed on the contralateral side of the articulating rod and securely tightened. The articulating rod was then removed and replaced with a second rod, augmented with a cross-connector over the T11 to T12 space for added stability. A chest tube was placed, and appropriate postoperative imaging was obtained to evaluate for pneumothorax. The patient was then closed in a standard fashion. IONM remained stable throughout the case. The total operative time was 478 minutes, with an estimated blood loss of 750 cc.

#### Postoperative Outcome

Thoracolumbar kyphotic angle was 35°, showing a 40° improvement from preoperative measurements. Anteroposterior and lateral intraoperative films showed adequate placement of hardware (Figure 1B and C). The patient was observed in our standard step-down unit until postoperative day 10, after which she was transferred to inpatient rehabilitation. No postoperative medical complications were noted, and our patient did not require a blood transfusion. Of note, due to the patient's significant preoperative opioid necessity, our acute pain management service was consulted for the comanagement of postoperative pain.

#### Follow-Up

At a 2-year postoperative evaluation, the patient had significant improvement in back and lower extremity pain, sagittal deformity, and claudication symptoms. No hardware issues were noted on follow-up imaging, and preservation of immediate postoperative correction was maintained (Figure 1D).

### DISCUSSION

Advancements in surgical technique and spinal technologies have allowed for high degrees of correction in sagittal deformity of various etiologies. In our patient who presented with delayed posttraumatic TLK, her severe (>70°) focal kyphosis necessitated extreme focal correction to allow for the return of sagittal alignment. The greatest concern for any deformity correction, especially in cases of hyperkyphosis of the distal thoracic spine, includes injury to neural elements due to cord distraction/ compression during the corrective maneuvers.<sup>6</sup> In addition, we believed that the use of pedicle subtraction osteotomy (PSO), while they may have provided the necessary degrees of correction, would have increased the overall risk profile with minimal gain considering her deformity was not rigid.<sup>7</sup>

Therefore, we opted to provide anterior support with a lateral interbody device, maximize lordotic potential with an anterior column release, and minimize the osteotomy necessary to achieve the necessary correction. However, the concern still remained of neurological compromise during the maneuver to distract the disc space open to allow for the interbody device to be introduced in a highly collapsed space and the maneuver to compress the thoracolumbar spine into the proper alignment. Here, the articulating rod was used in conjunction with IONM to safely achieve the planned correction. Depending on the points of fixation, we are able to directly control compression and distraction of the spine with rod holders on the articulating rod without the need for in situ benders. Additionally, any IONM changes can be detected instantaneously and can thereby be reversed quickly.




Figure 3. Illustration of placed articulating rod. (A) With the articulating rod in place within the pedicle screw tulip and all set caps finally tightened, downward pressure utilizing heavy rod holders allows for distraction of the disc space immediately inferior to the central pivot point. This allows for increased space to introduce an expandable lateral interbody device safely. (B) The cranial set caps are then loosened. When downward pressure is applied now, it allows for compression. With the sequential expansion of the interbody device with a continued downward force vector, the construct and spine can be placed in mild extension. Intraoperative neuromonitoring during these maneuvers is critical.

The intraoperative video shows the active use of the articulating rod (Supplemental Video 1). By placing the articulating rod with the central pivot hinge over the kyphotic point and finally tightening it with set caps, we can distract the disc space safely (Figure 3A). This allowed us to safely introduce our lateral interbody device (Supplemental Video 2). When we fix the articulating rod to the pedicle screws at the caudal end of the construct without final tightening set caps at the proximal end, we are able to compress the spine. With downward pressure, the articulating rod slides within the tulip heads (Figure 3B). We combined this maneuver with the expansion of the interbody device to achieve an adequate correction—approximately 40° comparable to a standard PSO.<sup>8</sup>

#### Limitations

While this is the first reported use of an articulating rod in literature for complex thoracolumbar correction of sagittal deformity, the applicability of this technology will need to be determined with additional trials. However, its inaugural use without complication, even at 2-year follow-up, is encouraging. Future reports should include long-term functional outcomes to provide a more complete patient assessment.

## CONCLUSION

Our use of a custom-made articulating rod allowed for the safe placement of an expandable lateral interbody device and controlled compression and distraction during a severe TLK correction case. We were able to achieve a 40° correction without the need for multilevel PSOs or a vertebral column resection, which we believe contributed to reduced risk and ultimate operative success. While more trials are necessary prior to further implementation, the use of articulating rods in select patients requiring extensive kyphotic corrections is technically feasible and may offer benefits in risk reduction during corrective maneuvers.

### REFERENCES

1. Chen J, Shao X-X, Sui W-Y, et al. Risk factors for neurological complications in severe and rigid spinal deformity correction of 177 cases. *BMC Neurol*. 2020;20(1):433. doi:10.1186/s12883-020-02012-8

2. Smith JS, Shaffrey CI, Glassman SD, et al. Risk-benefit assessment of surgery for adult scoliosis: an analysis based on patient age. *Spine (Phila Pa 1986)*. 1976;36(10):817–824. doi:10.1097/BRS.0b013e3181e21783

3. Wang X, Lenke LG, Thuet E, Blanke K, Koester LA, Roth M. Deformity angular ratio describes the severity of spinal deformity and predicts the risk of neurologic deficit in posterior vertebral column resection surgery. *Spine (Phila Pa 1986)*. 2016;41(18):1447–1455. doi:10.1097/BRS.000000000001547

4. Kelly MP, Lenke LG, Shaffrey CI, et al. Evaluation of complications and neurological deficits with three-column spine reconstructions for complex spinal deformity: a retrospective scoli-RISK-1 study. *Neurosurg Focus*. 2014;36(5):E17. doi:10.3171/2014.2.FO-CUS1419

5. Schwab FJ, Hawkinson N, Lafage V, et al. Risk factors for major peri-operative complications in adult spinal deformity surgery: a multi-center review of 953 consecutive patients. *Eur Spine J*. 2012;21(12):2603–2610. doi:10.1007/s00586-012-2370-4

6. Biscevic M, Sehic A, Krupic F. Intraoperative neuromonitoring in spine deformity surgery: modalities, advantages, limitations, medicolegal issues - surgeons' views. *EFORT Open Rev.* 2020;5(1):9–16. doi:10.1302/2058-5241.5.180032

7. Joiner EF, Mummaneni PV, Shaffrey CI, Chan AK. Posteriorbased osteotomies for deformity correction. *Neurosurg Clin N Am.* 2023;34(4):555–566. doi:10.1016/j.nec.2023.06.002

8. Berven SH, Deviren V, Smith JA, et al. Management of fixed sagittal plane deformity: results of the transpedicular wedge resection osteotomy. *Spine (Phila Pa 1986)*. 1976;26(18):2036–2043. doi:10.1097/00007632-200109150-00020

**Funding:** The authors received no financial support for the research, authorship, and/or publication of this article.

# **Declaration of Conflicting Interests:** The authors report no conflicts of interest in this work.

**Disclosures:** Timur Urakov, MD, discloses that he is a consultant for Globus Medical, Inc. The remaining authors have no other personal, financial, or institutional interests in any of the drugs, materials, or devices described in this article.

**Ethics Statement:** No Institutional Review Board approval was necessary per institutional policy. Consent for use of de-identified clinical information for publication was obtained.

**Corresponding Author:** Timur Urakov, Department of Neurological Surgery, University of Miami Miller School of Medicine, 1095 NW 14th Terrace, Miami, FL 33136, USA; turakov@med. miami.edu

Copyright © 2025 ISASS. The IJSS is an open access journal following the Creative Commons Licensing Agreement CC BY-NC-ND. To learn more or order reprints, visit http:// ijssurgery.com.